

For Research Use Only

MMP14 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67225-3-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 67225-3-PBS | GenBank Accession Number: BC064803 | Purification Method: Protein G purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 4323 | CloneNo.: 2D3B8 |
| Source: Mouse | UNIPROT ID: P50281 | |
| Isotype: IgG1 | Full Name: matrix metalloproteinase 14 (membrane-inserted) | |
| Immunogen Catalog Number: AG6062 | Calculated MW: 66 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

67225-3-PBS targets MMP14 as part of a matched antibody pair:

MP51349-1: 67225-2-PBS capture and 67225-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51349-1, MMP14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67225-2-PBS. Detection antibody: 67225-3-PBS.



Cytometric bead array standard curve of MP51349-1, MMP14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67225-2-PBS. Detection antibody: 67225-3-PBS. Standard:Ag6062. Range: 0.098-100 ng/mL.