

For Research Use Only

# ABTB1 Monoclonal antibody

Catalog Number: 67246-1-Ig



## Basic Information

Catalog Number:

67246-1-Ig

Size:

150ul , Concentration: 1000 µg/ml by  
Bradford method using BSA as the  
standard;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28959

GenBank Accession Number:

BC011858

GeneID (NCBI):

80325

UNIPROT ID:

Q969K4

Full Name:

ankyrin repeat and BTB (POZ) domain  
containing 1

Calculated MW:

54 kDa

Purification Method:

Protein G purification

CloneNo.:

2D2B11

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, pig

Positive Controls:

WB : pig cerebellum tissue,

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



pig cerebellum tissue were subjected to SDS PAGE followed by western blot with 67246-1-Ig (ABTB1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.