

For Research Use Only

ABTB1 Monoclonal antibody

Catalog Number: 67246-1-Ig



Basic Information

Catalog Number: 67246-1-Ig	GenBank Accession Number: BC011858	Purification Method: Protein G purification
Size: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 80325	CloneNo.: 2D2B11
Source: Mouse	Full Name: ankyrin repeat and BTB (POZ) domain containing 1	Recommended Dilutions:
Isotype: IgG1	Calculated MW: 54 kDa	
Immunogen Catalog Number: AG28959		

Applications

Tested Applications: WB,ELISA	Positive Controls:
Species Specificity: Human, pig	WB : pig cerebellum tissue,

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig cerebellum tissue were subjected to SDS PAGE followed by western blot with 67246-1-Ig (ABTB1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.