

For Research Use Only

# CD5 Monoclonal antibody, PBS Only

Catalog Number: 67278-1-PBS



## Basic Information

<b>Catalog Number:</b> 67278-1-PBS	<b>GenBank Accession Number:</b> BC027901	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 921	<b>CloneNo.:</b> 4D12C10
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000110448	
<b>Isotype:</b> IgG2b	<b>UNIPROT ID:</b> P06127	
<b>Immunogen Catalog Number:</b> AG11065	<b>Full Name:</b> CD5 molecule	
	<b>Calculated MW:</b> 495 aa, 55 kDa	

## Applications

**Tested Applications:**  
IHC, IF-P, Indirect ELISA

**Species Specificity:**  
human

## Background Information

CD5 is a type I transmembrane glycoprotein of the scavenger receptor cysteine-rich family (PMID: 12403363). CD5 is expressed on a majority of thymocytes, mature T cells, B cell subsets, and peripheral blood dendritic cells (PMID: 9858516; 6156984; 10379049; 1384337). CD5 may act as a receptor in regulating T-cell proliferation. It functions as a negative regulator of TCR signaling during thymocyte development (PMID: 7542801).

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67278-1-Ig (CD5 antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67278-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67278-1-Ig (CD5 antibody) at dilution of 1:2000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67278-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CD5 antibody (67278-1-Ig, Clone: 4D12C10) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67278-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CD5 antibody (67278-1-Ig, Clone: 4D12C10) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67278-1-PBS in a different storage buffer formulation.