

For Research Use Only

DLX3 Monoclonal antibody

Catalog Number: 67282-1-Ig



Basic Information

Catalog Number:

67282-1-Ig

GenBank Accession Number:

BC028970

Purification Method:

Protein A purification

Size:

150ul, Concentration: 1900 ug/ml by
Nanodrop and 1000 ug/ml by Bradford
method using BSA as the standard;

GeneID (NCBI):

1747

CloneNo.:

1H8E1

UNIPROT ID:

O60479

Recommended Dilutions:

WB 1:2000-1:10000

Source:

Mouse

Full Name:

distal-less homeobox 3

Isotype:

IgG2a

Calculated MW:

287 aa, 32 kDa

Immunogen Catalog Number:

AG29313

Observed MW:

50 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : pig brain tissue, HEK-293 cells, JAR cells, NCCIT
cells, rat brain tissue, mouse brain tissue, mouse
cerebellum tissue

Species Specificity:

Human, mouse, rat, pig

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67282-1-Ig (DLX3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.