

For Research Use Only

ALG13 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67283-1-PBS



Basic Information

Catalog Number:

67283-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14955

GenBank Accession Number:

BC005336

GeneID (NCBI):

79868

UNIPROT ID:

Q9NP73

Full Name:

asparagine-linked glycosylation 13 homolog

Calculated MW:

165 aa, 18 kDa

Observed MW:

130-135 kDa

Purification Method:

Protein G purification

CloneNo.:

3B8G3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67283-1-Ig (ALG13 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67283-1-PBS in a different storage buffer formulation.