

For Research Use Only

BRCA2 Monoclonal antibody

Catalog Number: 67285-1-Ig 1 Publications



Basic Information

Catalog Number: 67285-1-Ig	GenBank Accession Number: NM_000059	Purification Method: Protein G purification
Size: 150ul , Concentration: 553 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 675	CloneNo.: 3G1C7
Source: Mouse	UNIPROT ID: P51587	
Isotype: IgG1	Full Name: breast cancer 2, early onset	
	Calculated MW: 384 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Notable Publications

Author	Pubmed ID	Journal	Application
Ishraq Alim	35616648	Genet Med	

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

