

For Research Use Only

DAAM1 Monoclonal antibody, PBS Only

Catalog Number: 67287-1-PBS



Basic Information

Catalog Number:

67287-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6809

GenBank Accession Number:

BC064999

GeneID (NCBI):

23002

UNIPROT ID:

Q9Y4D1

Full Name:

dishevelled associated activator of morphogenesis 1

Calculated MW:

123 kDa

Observed MW:

120 kDa

Purification Method:

Protein A purification

CloneNo.:

3F11F11

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Dishevelled-associated activator of morphogenesis 1 (Daam1) is a member of the formin protein family that is implicated in actin assembly. Daam1 also plays a critical role in mediating Wnt/Fz signaling pathway and regulates cell polarity during development and tissue homeostasis. This antibody well recognizes the endogenous DAAM1 protein around 120 kDa.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67287-1-Ig (DAAM1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67287-1-PBS in a different storage buffer formulation.

