

For Research Use Only

ARL2BP Monoclonal antibody

Catalog Number: 67316-1-Ig



Basic Information

Catalog Number:

67316-1-Ig

GenBank Accession Number:

BC003087

Purification Method:

Protein G purification

Size:

150ul, Concentration: 1600 ug/ml by Nanodrop and 1000 ug/ml by Bradford method using BSA as the standard;

GeneID (NCBI):

23568

CloneNo.:

1G5C2

UNIPROT ID:

Q9Y2Y0

Recommended Dilutions:

WB 1:1000-1:6000

Source:

Mouse

Full Name:

ADP-ribosylation factor-like 2 binding protein

Isotype:

IgG1

Calculated MW:

19 kDa

Immunogen Catalog Number:

AG29261

Observed MW:

21 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : A549 cells, HeLa cells, Jurkat cells, HEK-293 cells, HepG2 cells

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.