

For Research Use Only

ARL2BP Monoclonal antibody, PBS Only



Catalog Number: 67316-1-PBS

Basic Information

Catalog Number: 67316-1-PBS	GenBank Accession Number: BC003087	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 23568	CloneNo.: 1G5C2
Source: Mouse	UNIPROT ID: Q9Y2Y0	
Isotype: IgG1	Full Name: ADP-ribosylation factor-like 2 binding protein	
Immunogen Catalog Number: AG29261	Calculated MW: 19 kDa	
	Observed MW: 21 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67316-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 67316-1-Ig (ARL2BP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67316-1-PBS in a different storage buffer formulation.