

For Research Use Only

TFE3 Monoclonal antibody, PBS Only

Catalog Number: 67319-1-PBS



Basic Information

Catalog Number:

67319-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6320

GenBank Accession Number:

BC026027

GeneID (NCBI):

7030

UNIPROT ID:

P19532

Full Name:

transcription factor binding to IGHM enhancer 3

Calculated MW:

62 kDa

Observed MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

1E3D7

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

Human, mouse, rat, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67319-1-Ig (TFE3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67319-1-PBS in a different storage buffer formulation.

