

For Research Use Only

NICN1 Monoclonal antibody, PBS Only



Catalog Number: 67328-1-PBS

Basic Information

Catalog Number:

67328-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG6581

GenBank Accession Number:

BC050653

GeneID (NCBI):

84276

UNIPROT ID:

Q9BSH3

Full Name:

nicolin 1

Calculated MW:

24 kDa

Observed MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

5B7G6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



U-251 cells were subjected to SDS PAGE followed by western blot with 67328-1-Ig (NICN1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67328-1-PBS in a different storage buffer formulation.