

For Research Use Only

# N-PAC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67340-3-PBS



## Basic Information

<b>Catalog Number:</b> 67340-3-PBS	<b>GenBank Accession Number:</b> BC064940	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 84656	<b>CloneNo.:</b> 2B1D2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q49A26	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> cytokine-like nuclear factor n-pac	
<b>Immunogen Catalog Number:</b> AG7133	<b>Calculated MW:</b> 61 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

67340-3-PBS targets N-PAC as part of a matched antibody pair.

MP51700-1: 67340-2-PBS capture and 67340-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51700-1, N-PAC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67340-2-PBS. Detection antibody: 67340-3-PBS.



Cytometric bead array standard curve of MP51700-1, N-PAC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67340-2-PBS. Detection antibody: 67340-3-PBS. Standard: Ag7133. Range: 1.563-200 ng/mL.