

For Research Use Only

ARL4C Monoclonal antibody

Catalog Number: 67356-1-Ig



Basic Information

Catalog Number: 67356-1-Ig	GenBank Accession Number: BC001051	Purification Method: Protein G purification
Size: 150ul	GeneID (NCBI): 10123	CloneNo.: 1B8B5
Source: Mouse	Full Name: ADP-ribosylation factor-like 4C	
Isotype: IgG1	Calculated MW: 21 kDa	
Immunogen Catalog Number: AG29496		

Applications

Tested Applications:
ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data