

For Research Use Only

BRD4 Monoclonal antibody

Catalog Number: 67374-1-Ig **2 Publications**



Basic Information

Catalog Number: 67374-1-Ig	GenBank Accession Number: BC035266	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23476	CloneNo.: 1F11A6
Source: Mouse	UNIPROT ID: O60885	
Isotype: IgG1	Full Name: bromodomain containing 4	
Immunogen Catalog Number: AG27545	Calculated MW: 1362 aa, 152 kDa	
	Observed MW: 150 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IF

Species Specificity:
Human

Cited Species:
human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Yanwen Lan	32596881	J Cell Mol Med	WB,IF
Zhixin Qiu	34615656	Cancer Discov	WB

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

