

For Research Use Only

RBMS2 Monoclonal antibody, PBS Only



Catalog Number: 67395-1-PBS

Basic Information

Catalog Number:

67395-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27721

GenBank Accession Number:

BC027863

GeneID (NCBI):

5939

UNIPROT ID:

Q15434

Full Name:

RNA binding motif, single stranded interacting protein 2

Calculated MW:

407 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

1E5C5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67395-1-Ig (RBMS2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67395-1-PBS in a different storage buffer formulation.