

For Research Use Only

GATAD2B Monoclonal antibody

Catalog Number: 67396-1-Ig



Basic Information

Catalog Number:

67396-1-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22615

GenBank Accession Number:

BC069419

GeneID (NCBI):

57459

UNIPROT ID:

Q8WXL9

Full Name:

GATA zinc finger domain containing
2B

Purification Method:

Protein G purification

CloneNo.:

2C3F5

Applications

Tested Applications:

ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

