

For Research Use Only

Secretogranin V Monoclonal antibody, PBS Only



Catalog Number: 67406-1-PBS

Basic Information

Catalog Number:

67406-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28607

GenBank Accession Number:

BC005349

GeneID (NCBI):

6447

UNIPROT ID:

P05408

Full Name:

secretogranin V (7B2 protein)

Calculated MW:

24 kDa

Purification Method:

Protein G purification

CloneNo.:

1B9H6

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pituitary adenoma tissue slide using 67406-1-Ig (Secretogranin V antibody) at dilution of 1:4000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67406-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human pituitary adenoma tissue slide using 67406-1-Ig (Secretogranin V antibody) at dilution of 1:4000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67406-1-PBS in a different storage buffer formulation.