

For Research Use Only

ERO1L Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67416-2-PBS



Basic Information

Catalog Number: 67416-2-PBS	GenBank Accession Number: BC008674	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 30001	CloneNo.: 1G10C6
Source: Mouse	UNIPROT ID: Q96HE7	
Isotype: IgG1	Full Name: ERO1-like (S. cerevisiae)	
Immunogen Catalog Number: AG29910	Calculated MW: 468 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67416-2-PBS targets ERO1L as part of a matched antibody pair.

MP51017-1: 67416-2-PBS capture and 67416-3-PBS detection (validated in Cytometric bead array)

MP51017-2: 67416-2-PBS capture and 67416-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51017-1, ERO1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67416-2-PBS. Detection antibody: 67416-3-PBS. Standard:Ag29910. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51017-2, ERO1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67416-2-PBS. Detection antibody: 67416-4-PBS. Standard:Ag29910. Range: 0.195-100 ng/mL