

For Research Use Only

AHCYL2 Monoclonal antibody, PBS Only



Catalog Number: 67484-1-PBS

Basic Information

Catalog Number:

67484-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19614

GenBank Accession Number:

BC024325

GeneID (NCBI):

23382

UNIPROT ID:

Q96HN2

Full Name:

S-adenosylhomocysteine hydrolase-like 2

Calculated MW:

611 aa, 67 kDa

Observed MW:

67 and 57 kDa

Purification Method:

Protein A purification

CloneNo.:

1D1F6

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 67484-1-Ig (AHCYL2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67484-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 67484-1-Ig (AHCYL2 antibody) at dilution of 1:600 (under 10x lens).. This data was developed using the same antibody clone with 67484-1-PBS in a different storage buffer formulation.