

For Research Use Only

RHOBTB3 Monoclonal antibody, PBS Only

Catalog Number: 67502-1-PBS



Basic Information

Catalog Number:

67502-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5782

GenBank Accession Number:

BC041337

GeneID (NCBI):

22836

UNIPROT ID:

O94955

Full Name:

Rho-related BTB domain containing 3

Calculated MW:

611 aa, 69 kDa

Observed MW:

69 kDa

Purification Method:

Protein G purification

CloneNo.:

1H2E3

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 67502-1-Ig (RHOBTB3 antibody) at dilution of 1:300 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67502-1-PBS in a different storage buffer formulation.



K-562 cells were subjected to SDS PAGE followed by western blot with 67502-1-Ig (RHOBTB3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67502-1-PBS in a different storage buffer formulation.