

For Research Use Only

RHOBTB3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67502-5-PBS



Basic Information

Catalog Number: 67502-5-PBS	GenBank Accession Number: BC041337	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 22836	CloneNo.: 2B12A6
Source: Mouse	UNIPROT ID: O94955	
Isotype: IgG1	Full Name: Rho-related BTB domain containing 3	
Immunogen Catalog Number: AG5782	Calculated MW: 611 aa, 69 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67502-5-PBS targets RHOBTB3 as part of a matched antibody pair:

MP51247-3: 67502-5-PBS capture and 67502-6-PBS detection (validated in Cytometric bead array)

MP51247-4: 67502-5-PBS capture and 67502-7-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51247-4, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-5-PBS. Detection antibody: 67502-7-PBS. Standard:Ag5782. Range: 1.563-100 ng/mL

Cytometric bead array standard curve of MP51247-3, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-5-PBS. Detection antibody: 67502-6-PBS. Standard:Ag5782. Range: 3.125-100 ng/mL