

For Research Use Only

# RBM8A,Y14 Monoclonal antibody, PBS Only (Detector)

Catalog Number:67541-3-PBS



## Basic Information

Catalog Number:

67541-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6800

GenBank Accession Number:

BC017088

GeneID (NCBI):

9939

UNIPROT ID:

Q9Y5S9

Full Name:

RNA binding motif protein 8A

Calculated MW:

20 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A5B6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Product Information

67541-3-PBS targets RBM8A,Y14 as part of a matched antibody pair:

MP51528-1: 67541-2-PBS capture and 67541-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

