

For Research Use Only

BTN3A3 Monoclonal antibody, PBS Only



Catalog Number: 67560-1-PBS

Basic Information

Catalog Number: 67560-1-PBS	GenBank Accession Number: BC015815	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10384	CloneNo.: 1A3B5
Source: Mouse	UNIPROT ID: O00478	
Isotype: IgG1	Full Name: butyrophilin, subfamily 3, member A3	
Immunogen Catalog Number: AG12755	Calculated MW: 584 aa, 65 kDa	
	Observed MW: 65 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67560-1-Ig (BTN3A3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67560-1-PBS in a different storage buffer formulation.