

For Research Use Only

TAF6L Monoclonal antibody, PBS Only



Catalog Number: 67569-1-PBS

Basic Information

Catalog Number:

67569-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG8202

GenBank Accession Number:

BC008785

GeneID (NCBI):

10629

UNIPROT ID:

Q9Y6J9

Full Name:

TAF6-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa

Calculated MW:

622 aa, 68 kDa

Observed MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

1G5C4

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67569-1-Ig (TAF6L antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67569-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 67569-1-Ig (TAF6L antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67569-1-PBS in a different storage buffer formulation.