

For Research Use Only

Vaspin Monoclonal antibody, PBS Only

Catalog Number: 67573-1-PBS



Basic Information

Catalog Number:

67573-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30136

GenBank Accession Number:

BC040857

GeneID (NCBI):

145264

UNIPROT ID:

Q8IW75

Full Name:

serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 12

Calculated MW:

414 aa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Protein A purification

CloneNo.:

2A1B6

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human adipose tissue were subjected to SDS PAGE followed by western blot with 67573-1-Ig (Vaspin antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67573-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Vaspin antibody (67573-1-Ig, Clone: 2A1B6) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 67573-1-PBS in a different storage buffer formulation.