

For Research Use Only

Galectin 10 Monoclonal antibody, PBS Only



Catalog Number: 67581-1-PBS

Basic Information

Catalog Number: 67581-1-PBS	GenBank Accession Number: BC119711	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1178	CloneNo.: 5D9B7
Source: Mouse	UNIPROT ID: Q05315	
Isotype: IgG1	Full Name: Charcot-Leyden crystal protein	
Immunogen Catalog Number: AG18569	Calculated MW: 142 aa, 16 kDa	
	Observed MW: 16 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HL-60 cells were subjected to SDS PAGE followed by western blot with 67581-1-Ig (CLC antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67581-1-PBS in a different storage buffer formulation.