

For Research Use Only

CALU Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67585-2-PBS



Basic Information

Catalog Number: 67585-2-PBS	GenBank Accession Number: BC013383	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 813	CloneNo.: 1C1D4
Source: Mouse	UNIPROT ID: O43852	
Isotype: IgG1	Full Name: calumenin	
Immunogen Catalog Number: AG12073	Calculated MW: 315 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67585-2-PBS targets CALU as part of a matched antibody pair:

MP50698-1: 67585-2-PBS capture and 67585-3-PBS detection (validated in Cytometric bead array)

MP50698-2: 67585-2-PBS capture and 67585-1-PBS detection (validated in Cytometric bead array)

MP50698-3: 67585-2-PBS capture and 67585-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50698-1, CALU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67585-2-PBS. Detection antibody: 67585-3-PBS. Standard:Ag12073. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50698-2, CALU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67585-2-PBS. Detection antibody: 67585-1-PBS. Standard:Ag12073. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50698-3, CALU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67585-2-PBS. Detection antibody: 67585-4-PBS. Standard:Ag12073. Range: 0.195-100 ng/mL