

For Research Use Only

Mitoferrin 1 Monoclonal antibody, PBS Only



Catalog Number: 67593-1-PBS

Basic Information

Catalog Number: 67593-1-PBS	GenBank Accession Number: BC132799	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51312	CloneNo.: 1B1B2
Source: Mouse	UNIPROT ID: Q9NYZ2	
Isotype: IgG1	Full Name: solute carrier family 25, member 37	
Immunogen Catalog Number: AG26216	Observed MW: 40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



TF-1 cells were subjected to SDS PAGE followed by western blot with 67593-1-Ig (Mitoferrin 1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67593-1-PBS in a different storage buffer formulation.



human placenta tissue were subjected to SDS PAGE followed by western blot with 67593-1-Ig (Mitoferrin 1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67593-1-PBS in a different storage buffer formulation.