

For Research Use Only

# CD226 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67653-2-PBS

## Basic Information

<b>Catalog Number:</b> 67653-2-PBS	<b>GenBank Accession Number:</b> BC074787	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100ug , Concentration: 1mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 10666	<b>CloneNo.:</b> 2D1A4
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000150637	
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> Q15762	
<b>Immunogen Catalog Number:</b> AG12403	<b>Full Name:</b> CD226 molecule	
	<b>Calculated MW:</b> 336 aa, 39 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Product Information

67653-2-PBS targets CD226 as part of a matched antibody pair:

MP50115-1: 67653-2-PBS capture and 67653-3-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50115-1, CD226 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67653-2-PBS. Detection antibody: 67653-3-PBS. Standard:Ag12403. Range: 0.391-100 ng/mL.