

For Research Use Only

LAMB1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67705-3-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 67705-3-PBS | GenBank Accession Number: BC113455 | Purification Method: Protein G Magarose purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 3912 | CloneNo.: 2C1B4 |
| Source: Mouse | UNIPROT ID: P07942 | |
| Isotype: IgG1 | Full Name: laminin, beta 1 | |
| Immunogen Catalog Number: AG19182 | Calculated MW: 1786 aa, 198 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67705-3-PBS targets LAMB1 as part of a matched antibody pair:

MP51088-1: 67705-2-PBS capture and 67705-3-PBS detection (validated in Cytometric bead array)

MP51088-2: 67705-4-PBS capture and 67705-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.


Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP51088-1, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-2-PBS. Detection antibody: 67705-3-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP51088-1, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-2-PBS. Detection antibody: 67705-3-PBS. Standard:Ag19182. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51088-2, LAMB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67705-4-PBS. Detection antibody: 67705-3-PBS. Standard:Ag19182. Range: 0.195-50 ng/mL