

For Research Use Only

GIMAP2 Monoclonal antibody

Catalog Number: 67720-1-Ig



Basic Information

Catalog Number: 67720-1-Ig	GenBank Accession Number: BC032345	Purification Method: Protein G purification
Size: 150ul , Concentration: 1500 µg/ml by Nanodrop and 749 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 26157	CloneNo.: 1F2A2
Source: Mouse	Full Name: GTPase, IMAP family member 2	Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG1	Calculated MW: 337 aa, 38 kDa	
Immunogen Catalog Number: AG4368	Observed MW: 33-38 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MOLT-4 cells, Jurkat cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MOLT-4 cells were subjected to SDS PAGE followed by western blot with 67720-1-Ig (GIMAP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.