

For Research Use Only

GIMAP2 Monoclonal antibody

Catalog Number: 67720-1-Ig



Basic Information

Catalog Number:

67720-1-Ig

Size:

150ul , Concentration: 1500 ug/ml by Nanodrop and 749 ug/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4368

GenBank Accession Number:

BC032345

GeneID (NCBI):

26157

UNIPROT ID:

Q9UG22

Full Name:

GTPase, IMAP family member 2

Calculated MW:

337 aa, 38 kDa

Observed MW:

33-38 kDa

Purification Method:

Protein G purification

CloneNo.:

1F2A2

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : MOLT-4 cells, Jurkat cells

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MOLT-4 cells were subjected to SDS PAGE followed by western blot with 67720-1-Ig (GIMAP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.