

For Research Use Only

VASH2 Monoclonal antibody, PBS Only

Catalog Number: 67753-1-PBS



Basic Information

Catalog Number: 67753-1-PBS	GenBank Accession Number: BC028194	Purification Method: Protein G purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 79805	CloneNo.: 3A11D10
Source: Mouse	UNIPROT ID: Q86V25	
Isotype: IgG1	Full Name: vasohibin 2	
Immunogen Catalog Number: AG29985	Calculated MW: 355 aa, 40 kDa	
	Observed MW: 30-34 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

VASH2, also named as vasohibin 2, is an important pro-angiogenesis factor in solid tumor. It has been reported that VASH2 is expressed in mononuclear cells mobilized from bone marrow to promote angiogenesis. VASH2 also plays a key role in axon formation. VASH2 is localized as a nuclear type (with 311 amino acid residues) and cytoplasmic type (with 355 amino acid residues and low abundance). Cytoplasmic VASH2 is associated with carcinoma angiogenesis, while nuclear VASH2 may be associated with cell proliferation. 67753-1-Ig antibody recognizes both the nuclear and cytoplasmic isoforms. (PMID:23615928; 19204325; 31235911; 26177649)

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67753-1-Ig (VASH2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 67753-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 67753-1-Ig (VASH2 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67753-1-PBS in a different storage buffer formulation.