

For Research Use Only

TMEM74 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67773-3-PBS



Basic Information

Catalog Number: 67773-3-PBS	GenBank Accession Number: BC030710	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 157753	CloneNo.: 2F3B12
Source: Mouse	UNIPROT ID: Q96NL1	
Isotype: IgG1	Full Name: transmembrane protein 74	
Immunogen Catalog Number: AG22359	Calculated MW: 305 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67773-3-PBS targets TMEM74 as part of a matched antibody pair:

MP51234-1: 67773-2-PBS capture and 67773-3-PBS detection (validated in Cytometric bead array)

MP51234-2: 67773-3-PBS capture and 67773-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.


Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP51234-1, TMEM74 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67773-2-PBS. Detection antibody: 67773-3-PBS. Standard:Ag22359. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51234-2, TMEM74 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67773-3-PBS. Detection antibody: 67773-4-PBS. Standard:Ag22359. Range: 0.391-100 ng/mL