

For Research Use Only

RAB11B Monoclonal antibody, PBS Only



Catalog Number: 67780-1-PBS

Basic Information

Catalog Number: 67780-1-PBS	GenBank Accession Number: BC110081	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9230	CloneNo.: 1G7D6
Source: Mouse	UNIPROT ID: Q15907	
Isotype: IgG1	Full Name: RAB11B, member RAS oncogene family	
Immunogen Catalog Number: AG25429	Calculated MW: 24 kDa	
	Observed MW: 24 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67780-1-Ig (RAB11B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67780-1-PBS in a different storage buffer formulation.