

For Research Use Only

RP2 Monoclonal antibody, PBS Only

Catalog Number: 67795-1-PBS



Basic Information

Catalog Number:

67795-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5630

GenBank Accession Number:

BC043348

GeneID (NCBI):

6102

UNIPROT ID:

O75695

Full Name:

retinitis pigmentosa 2 (X-Linked recessive)

Calculated MW:

350 aa, 40 kDa

Observed MW:

37-40 kDa

Purification Method:

Protein G purification

CloneNo.:

1G2A9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, rabbit

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 67795-1-Ig (RP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67795-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 67795-1-Ig (RP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67795-1-PBS in a different storage buffer formulation.