

For Research Use Only

# PKMYT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67806-4-PBS



## Basic Information

<b>Catalog Number:</b> 67806-4-PBS	<b>GenBank Accession Number:</b> NM_004203	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 9088	<b>CloneNo.:</b> 1G11F7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q99640	
<b>Isotype:</b> IgG2b	<b>Full Name:</b> protein kinase, membrane associated tyrosine/threonine 1	
<b>Immunogen Catalog Number:</b> AG30849	<b>Calculated MW:</b> 55 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

67806-4-PBS targets PKMYT1 as part of a matched antibody pair.

MP50951-2: 67806-4-PBS capture and 67806-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50951-2, PKMYT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67806-4-PBS. Detection antibody: 67806-5-PBS. Standard:Ag30849. Range: 0.781-100 ng/mL.