

For Research Use Only

SLC25A42 Monoclonal antibody, PBS Only



Catalog Number: 67819-1-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 67819-1-PBS | GenBank Accession Number: BC045598 | Purification Method: Protein A purification |
| Size: 100ug, Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 284439 | CloneNo.: 2F7G10 |
| Source: Mouse | UNIPROT ID: Q86VD7 | |
| Isotype: IgG2b | Full Name: solute carrier family 25, member 42 | |
| Immunogen Catalog Number: AG26379 | Calculated MW: 318 aa, 35 kDa | |
| | Observed MW: 33-35 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 67819-1-Ig (SLC25A42 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67819-1-PBS in a different storage buffer formulation.