

For Research Use Only

INADL Monoclonal antibody

Catalog Number: 67833-1-Ig



Basic Information

Catalog Number:

67833-1-Ig

Size:

150ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15953

GenBank Accession Number:

BC142661

GeneID (NCBI):

10207

UNIPROT ID:

Q8NI35

Full Name:

InaD-like (Drosophila)

Calculated MW:

1801 aa, 196 kDa

Observed MW:

196 kDa

Purification Method:

Protein G purification

CloneNo.:

3C6A4

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Rat

Positive Controls:

WB : LNCaP cells, HeLa cells, Jurkat cells, HEK-293
cells, K 562 cells, PC-12 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67833-1-Ig (INADL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.