

For Research Use Only

INADL Monoclonal antibody, PBS Only



Catalog Number: 67833-1-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 67833-1-PBS | GenBank Accession Number: BC142661 | Purification Method: Protein G purification |
| Size: 100ug, Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 10207 | CloneNo.: 3C6A4 |
| Source: Mouse | UNIPROT ID: Q8NI35 | |
| Isotype: IgG1 | Full Name: InaD-like (Drosophila) | |
| Immunogen Catalog Number: AG15953 | Calculated MW: 1801 aa, 196 kDa | |
| | Observed MW: 196 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Rat

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67833-1-Ig (INADL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67833-1-PBS in a different storage buffer formulation.