

For Research Use Only

INADL Monoclonal antibody, PBS Only



Catalog Number: 67833-1-PBS

Basic Information

Catalog Number:

67833-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15953

GenBank Accession Number:

BC142661

GeneID (NCBI):

10207

UNIPROT ID:

Q8NI35

Full Name:

InaD-like (Drosophila)

Calculated MW:

1801 aa, 196 kDa

Observed MW:

196 kDa

Purification Method:

Protein G purification

CloneNo.:

3C6A4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67833-1-Ig (INADL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67833-1-PBS in a different storage buffer formulation.