

For Research Use Only

# CAPNS1 Monoclonal antibody, PBS Only



Catalog Number: 67838-1-PBS

## Basic Information

Catalog Number:

67838-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21885

GenBank Accession Number:

BC011903

GeneID (NCBI):

826

UNIPROT ID:

P04632

Full Name:

calpain, small subunit 1

Calculated MW:

28 kDa

Purification Method:

Protein G purification

CloneNo.:

2B6G6

## Applications

Tested Applications:

WB, IF, Indirect ELISA

Species Specificity:

Human, rat

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A375 cells were subjected to SDS PAGE followed by western blot with 67838-1-Ig (CAPNS1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67838-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using CAPNS1 antibody (67838-1-Ig, Clone: 2B6G6) at dilution of 1:500 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67838-1-PBS in a different storage buffer formulation.