

For Research Use Only

SH3BP5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67841-5-PBS



Basic Information

Catalog Number: 67841-5-PBS	GenBank Accession Number: BC010123	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9467	CloneNo.: 3D8G8
Source: Mouse	UNIPROT ID: O60239	
Isotype: IgG1	Full Name: SH3-domain binding protein 5 (BTK-associated)	
Immunogen Catalog Number: AG30252	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

67841-5-PBS targets SH3BP5 as part of a matched antibody pair:

MP50694-2: 67841-4-PBS capture and 67841-5-PBS detection (validated in Cytometric bead array)

MP50694-3: 67841-2-PBS capture and 67841-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50694-2, SH3BP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67841-4-PBS. Detection antibody: 67841-5-PBS. Standard:Ag30252. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50694-3, SH3BP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67841-2-PBS. Detection antibody: 67841-5-PBS. Standard:Ag30252. Range: 0.391-50 ng/mL.