

For Research Use Only

PGA5 Monoclonal antibody

Catalog Number: 67848-1-Ig



Basic Information

Catalog Number: 67848-1-Ig	GenBank Accession Number: BC029055	Purification Method: Protein G purification
Size: 150ul , Concentration: 2287 µg/ml by Nanodrop;	GeneID (NCBI): 5222	CloneNo.: 1C5E12
Source: Mouse	UNIPROT ID: P0DJ09	
Isotype: IgG1	Full Name: pepsinogen 5, group I (pepsinogen A)	
Immunogen Catalog Number: AG11657	Calculated MW: 388 aa, 42 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

