

For Research Use Only

SCAMP3 Monoclonal antibody, PBS Only



Catalog Number: 67932-1-PBS

Basic Information

Catalog Number: 67932-1-PBS	GenBank Accession Number: BC000161	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10067	CloneNo.: 1C1E6
Source: Mouse	UNIPROT ID: O14828	
Isotype: IgG1	Full Name: secretory carrier membrane protein 3	
Immunogen Catalog Number: AG26073	Calculated MW: 38 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications:
WB, IF, ELISA

Species Specificity:
Human, rat, mouse

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67932-1-Ig (SCAMP3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67932-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using SCAMP3 antibody (67932-1-Ig, Clone: 1C1E6) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67932-1-PBS in a different storage buffer formulation.