

For Research Use Only

PUF-A, KIAA0020 Monoclonal antibody, PBS Only



Catalog Number: 67941-1-PBS

Basic Information

Catalog Number:

67941-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31324

GenBank Accession Number:

BC016137

GeneID (NCBI):

9933

UNIPROT ID:

Q15397

Full Name:

KIAA0020

Calculated MW:

648 aa, 74 kDa

Observed MW:

68-70 kDa

Purification Method:

Protein G purification

CloneNo.:

1G4B7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67941-1-Ig (KIAA0020 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67941-1-PBS in a different storage buffer formulation.