

For Research Use Only

SLC7A5 Monoclonal antibody, PBS Only



Catalog Number: 67951-1-PBS

Basic Information

Catalog Number:

67951-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG20291

GenBank Accession Number:

BC039692

GeneID (NCBI):

8140

UNIPROT ID:

Q01650

Full Name:

solute carrier family 7 (cationic amino acid transporter, y+ system), member 5

Calculated MW:

507 aa, 55 kDa

Observed MW:

35-50 kDa

Purification Method:

Protein A purification

CloneNo.:

2G5H3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67951-1-Ig (SLC7A5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67951-1-PBS in a different storage buffer formulation.