

For Research Use Only

RON, MST1R Monoclonal antibody, PBS Only



Catalog Number: 67954-1-PBS

Basic Information

Catalog Number:

67954-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28289

GenBank Accession Number:

BQ936677 NM_002447

GeneID (NCBI):

4486

UNIPROT ID:

Q04912

Full Name:

macrophage stimulating 1 receptor (c-
met-related tyrosine kinase)

Calculated MW:

152 kDa

Observed MW:

150-185 kDa

Purification Method:

Protein G purification

CloneNo.:

2D8D1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 67954-1-Ig (RON, MST1R antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67954-1-PBS in a different storage buffer formulation.