

For Research Use Only

# Caspase 7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67956-3-PBS



## Basic Information

|                                                             |                                                                      |                                                         |
|-------------------------------------------------------------|----------------------------------------------------------------------|---------------------------------------------------------|
| <b>Catalog Number:</b><br>67956-3-PBS                       | <b>GenBank Accession Number:</b><br>BC015799                         | <b>Purification Method:</b><br>Euglobulin precipitation |
| <b>Size:</b><br>100ug , Concentration: 1 mg/ml by Nanodrop; | <b>GeneID (NCBI):</b><br>840                                         | <b>CloneNo.:</b><br>1C1E3                               |
| <b>Source:</b><br>Mouse                                     | <b>UNIPROT ID:</b><br>P55210                                         |                                                         |
| <b>Isotype:</b><br>IgM                                      | <b>Full Name:</b><br>caspase 7, apoptosis-related cysteine peptidase |                                                         |
| <b>Immunogen Catalog Number:</b><br>AG27601                 | <b>Calculated MW:</b><br>303 aa, 34 kDa                              |                                                         |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

67956-3-PBS targets Caspase 7 as part of a matched antibody pair:

MP51724-1: 67956-2-PBS capture and 67956-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51724-1, Caspase 7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67956-2-PBS. Detection antibody: 67956-3-PBS.



Cytometric bead array standard curve of MP51724-1, Caspase 7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67956-2-PBS. Detection antibody: 67956-3-PBS. Standard:Ag27601. Range: 1.563-200 ng/mL.